




1	<b>Parathormone-related Peptide (PTHrP)</b> Test Code : 8397002 - Laboratory Services			
	Methodology : RIA - Performed : Tue 08:00 - Reported : 1 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma, Frozen	Lavender top (EDTA) (Trasylol)	2mL	<0°C
<b>Additional Information</b> 				
Obtain EDTA-Trasylol (APROTININ) tube from laboratory. Sample should be centrifuged within 30 mins of collection and frozen within 2 hours of collection. After centrifugation, the plasma is separated and transferred to the laboratory. The plasma can be stored and shipped at -20°C or colder for up to 4 weeks.				

2	<b>Zinc (Zn), Erythrocyte</b> Test Code : 8463002 - Laboratory Services			
	Method : ICP-MS - Performed : Tue 08:00 Fri 08:00 - Reported : Fri 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Whole Blood (EDTA)	Lavender top (EDTA)	2mL	2-8°C

3	<b>Copper (Cu), erythrocyte</b> Test Code : 8252504 - Laboratory Services			
	Method : ICP-MS - Performed : Tue 08:00 Fri 08:00 - Reported : Fri 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Whole Blood (EDTA)	Lavender top (EDTA)	2mL	2-8°C

4	<b>Magnesium (Mg), Erythrocyte</b> Test Code : 8373506 - Laboratory Services			
	Method : ICP-MS - Performed : Tue 08:00 Fri 08:00 - Reported : Fri 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Whole Blood (EDTA)	Lavender top (EDTA)	2mL	2-8°C

5	Selenium, RBCs			
	Test Code : 8425505 - Laboratory Services			
	Method : ICP-MS - Performed : Tue 08:00 Fri 08:00 - Reported : Fri 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Whole Blood (EDTA)	Lavender top (EDTA)	2mL	2-8°C

6	Copper (Cu), Tissue			
	Test Code : 8252503 - Laboratory Services			
	Method : ICP-MS - Performed : Tue 08:00 - Reported : 2 weeks			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Tissue Sample	Call the Laboratory	0,25gr	<0°C