


# Laboratory Test Updates

(effective on August 12, 2024)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8818083	Polycystic Ovarian Syndrome (PCOS) Diagnostic Tests						X				
2	C8631605	Neuron Specific Enolase (NSE)			X							
3	C8631610	Neuron Specific Enolase (NSE), CSF			X							
4	C8925620	Sexually Transmitted Disease Panel - 13 Parameters				X			X			
5	C8727202	Giardia Antigen			X							

1	<b>Polycystic Ovarian Syndrome (PCOS) Diagnostic Tests</b> Test Code : <b>8818083</b> - Laboratory Services			
	Methodology : - Performed : <b>Mon 08:00 Thu 08:00</b> - Reported : <b>Tue 18:00 Fri 18:00</b>			
	Specimen Type	Container	Volume	Transport
	Serum	Separator Tube	5mL	2-8°C
<b>Additional Information</b> 				
Panel Contains; Follicle Stimulating Hormone (FSH), Luteinizing Hormone (LH), 17-Alpha-Hydroxyprogesterone (Direct), Dehydroepiandrosterone Sulphate (DHEA-S), Testosterone, Total, Sex Hormone Binding Globulin (SHBG), Anti-Mullerian Hormone (AMH/MIS), Thyroid Stimulating Hormone (TSH), Insulin.				

2	<b>Neuron Specific Enolase (NSE)</b> Test Code : <b>8631605</b> - Laboratory Services			
	Method : <b>ECLIA</b> - Performed : <b>Mon 08:00 Fri 08:00</b> - Reported : <b>Mon 18:00 Fri 18:00</b>			
	Sample Name	Container	Amount	Transfer
Serum	Separator Tube	2mL	2-8°C	

3	<b>Neuron Specific Enolase (NSE), CSF</b> Test Code : <b>8631610</b> - Laboratory Services			
	Methodology : <b>ECLIA</b> - Performed : <b>Mon 08:00 Fri 08:00</b> - Reported : <b>Mon 18:00 Fri 18:00</b>			
	Specimen Type	Container	Volume	Transport
CSF	Sterile CSF Tube	1mL	2-8°C	

<b>4</b>	<b>Sexually Transmitted Disease Panel - 13 Parameters</b> Test Code : <b>8925620</b> - Laboratory Services			
	Method : - Performed : <b>Mon 08:00 Wed 08:00</b> - Reported : <b>Tue 18:00 Thu 18:00</b>			
	<b>Sample Name</b>	<b>Container</b>	<b>Amount</b>	<b>Transfer</b>
	Serum and Urine	Separator Tube and Plastic Collection Cup	2mL	2-8°C
<b>Additional Information</b> <span style="float: right;">▼</span>				
<p>This panel includes; Chlamydia trachomatis, Mycoplasma genitalium, Neisseria gonorrhoeae, Trichomonas vaginalis, Gardnerella vaginalis, Ureaplasma urealyticum/parvum, Herpes Simplex Virus 1/2, Mycoplasma hominis, HIV Antibody (1-2) + p24 Antigen, Hepatitis B s Antigen (HBsAg), Hepatitis C Antibody (anti-HCV), Syphilis Profile. All sample types defined in the test panel are required.</p>				
<b>Serum and Urine</b> The first morning urine sample is taken.				

<b>5</b>	<b>Giardia Antigen</b> Test Code : <b>8727202</b> - Laboratory Services			
	Methodology : <b>IMMUNOCHROMATOGRAPHY</b> - Performed : <b>Mon-Fri 08:00 Sat 08:00</b> - Reported : <b>Mon-Fri 18:00 Mon 18:00</b>			
	<b>Specimen Type</b>	<b>Container</b>	<b>Volume</b>	<b>Transport</b>
	Stool	Plastic Collection Cup	1gr	2-8°C