

Laboratory Test Updates

(effective on July 13, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8818061	Anti CD20 Treatment Monitoring, Bone Marrow							X			
2	C8818060	Anti CD20 Treatment Monitoring, Whole Blood							X			
3	C8818062	Anti CD20 Treatment Monitoring, Biologic Fluid							X			
4	C8925686	Lymphocyte Subsets, Peripheral Blood							X			
5	C8925687	Lymphocyte Subsets, Bone Marrow							X			
6	C8925688	Lymphocyte Subsets, Biologic Fluid							X			
7	C8636107	PNH (Paroxysmal Nocturnal Hemoglobinuria) Panel							X			
8	C8818064	Primary Immunodeficiency Panel							X			
9	C8202406	Paliperidone									X	
10	C8616113	Fibrosis-4 (FIB-4) Score						X				

1

Anti CD20 Treatment Monitoring, Bone Marrow

Test Code : 8818061 - Laboratory Services

Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 4 days

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Bone Marrow	Lavender top (EDTA)	5mL	18-24°C

Additional Information

Please contact with the laboratory before collecting the sample. Send specimen in original tube. Do not aliquot.

2

Anti CD20 Treatment Monitoring, Whole Blood

Test Code : 8818060 - Laboratory Services

Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 4 days

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Whole Blood	Lavender top (EDTA)	5mL	18-24°C

Additional Information

Please contact with the laboratory before collecting the sample. Send specimen in original tube.

3

Anti CD20 Treatment Monitoring, Biologic Fluid

Test Code : 8818062 - Laboratory Services

Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 4 days

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Biologic Fluid	Lavender top (EDTA)	5mL	18-24°C

Additional Information

Please contact with the laboratory before collecting the sample. Send specimen in original tube. Do not aliquot.

4	Lymphocyte Subsets, Peripheral Blood Test Code : 8925686 - Laboratory Services			
	Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 3 days			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (EDTA)	Lavender top (EDTA)	5mL	18-24°C
Additional Information 				
Please contact with the laboratory before collecting the sample.Fresh sample is required.It is recommended that specimens arrive within 24 hours of draw.Test contains CD3,CD4, CD8, CD19, CD16/56.				
Whole Blood (EDTA) Unacceptable Conditions: Hemolysis and Clotted specimens				

5	Lymphocyte Subsets, Bone Marrow Test Code : 8925687 - Laboratory Services			
	Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 3 days			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Bone Marrow	Lavender top (EDTA)	5mL	18-24°C
Additional Information 				
Please contact with the laboratory before collecting the sample.Fresh sample is required.It is recommended that specimens arrive within 24 hours of draw.Test contains CD3,CD4, CD8, CD19, CD16/56.				
Bone Marrow Unacceptable Conditions: Hemolysis and Clotted specimens				

6	Lymphocyte Subsets, Biologic Fluid Test Code : 8925688 - Laboratory Services			
	Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 3 days			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Biologic Fluid	Lavender top (EDTA)	5mL	18-24°C
Additional Information 				
Please contact with the laboratory before collecting the sample.Fresh sample is required.It is recommended that specimens arrive within 24 hours of draw.Test contains CD3,CD4, CD8, CD19, CD16/56.				
Biologic Fluid Unacceptable Conditions: Hemolysis and Clotted specimens				

7

PNH (Paroxysmal Nocturnal Hemoglobinuria) Panel

Test Code : 8636107 - Laboratory Services

Methodology : FLOW CYTOMETRY - Performed : Mon 08:00 Tue 08:00 Wed 08:00 Thu 08:00 Fri 08:00 Sat 08:00 -
Reported : Wed 18:00 Thu 18:00 Fri 18:00 Mon 18:00 Tue 18:00 Tue 18:00

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Whole Blood (EDTA)	Lavender top (EDTA)	5mL	18-24°C

Additional Information

Please contact with the laboratory before collecting the sample. Fresh sample is required. It is recommended that specimens arrive within 24 hours of draw.

Whole Blood (EDTA)

Unacceptable Conditions: Hemolysis and Clotted specimens

8

Primary Immunodeficiency Panel Test Code : 8818064 - Laboratory Services			
Methodology : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 4 days			
Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Whole Blood	Lavender top (EDTA)		
Whole Blood (EDTA)	Lavender top (EDTA) (CBC)	5mL	18-24°C
Additional Information 			
Please contact with the laboratory before collecting the sample. Send specimen in original tube. Do not aliquot.			

9	Paliperidone Test Code : 8202406 - Laboratory Services		
	Methodology : LC-MS/MS - Performed : Tue 08:00 - Reported : 1 week		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Plain Tube	1mL
	Transport		2-8°C

10	Fibrosis-4 (FIB-4) Score		
	Test Code : 8616113 - Laboratory Services		
	Methodology : - Performed : Mon-Sat 08:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum and Whole Blood (EDTA)	Separator + Lavender top (EDTA)	2mL
			Transport
			2-8°C