

Laboratory Test Updates

(effective on October 14, 2024)

| | Test Code | Test Name | Name Change | Methodology | Performed/Reported Schedule | Specimen Requirements | Reference Interval | New Tests | Additional Information | Internalized Tests | Reactivated Tests | Inactive |
|----|-----------|--|-------------|-------------|-----------------------------|-----------------------|--------------------|-----------|------------------------|--------------------|-------------------|----------|
| 1 | C8780122 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Gastric Fluid, Fasting | | | | X | | | | | | |
| 2 | C8780123 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Sputum | | | | X | | | | | | |
| 3 | C8780125 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), CSF | | | | X | | | | | | |
| 4 | C8780126 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Bronchoalveolar Lavage | | | | X | | | | | | |
| 5 | C8780131 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Cystic Fluid | | | | X | | | | | | |
| 6 | C8779946 | Aspergillus DNA (PCR), Sputum | | | | X | | | | | | |
| 7 | C8779950 | Aspergillus DNA (PCR), Pus | | | | X | | | | | | |
| 8 | C8748101 | Candida DNA (PCR) | | | | X | | | | | | |
| 9 | C8748601 | Chlamydia pneumoniae DNA (PCR) | | | | X | | | | | | |
| 10 | C8754201 | Legionella pneumophila DNA (PCR) | | | | X | | | | | | |
| 11 | C8779914 | Nocardia DNA (PCR) | | | | X | | | | | | |
| 12 | C8779823 | Respiratory Syncytial Virus (RSV) RNA (PCR) | | | | X | | | | | | |
| 13 | C8745001 | Respiratory Panel, PCR | | | | X | | | | | | |
| 14 | C8780132 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Paraffin Section | | | | X | | | | | | |

| | | | | | | | | | | | | |
|----|----------|--------------------------------------|--|---|---|--|--|--|--|---|--|--|
| 15 | C8623530 | Autoimmune Encephalitis Panel, Serum | | X | X | | | | | X | | |
| 16 | C8623531 | Autoimmune Encephalitis Panel, CSF | | X | X | | | | | X | | |

| 1 | Tuberculosis (<i>Mycobacterium tuberculosis</i>) DNA, (PCR), Gastric Fluid, Fasting Test Code : 8780122 - Laboratory Services | | | |
|---|--|---------------------|--------|----------|
| | Method : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Sample Name | Container | Amount | Transfer |
| | Gastric Fluid, Fasting | Biologic Fluid Tube | 2mL | 2-8°C |

| 2 | Tuberculosis (<i>Mycobacterium tuberculosis</i>) DNA, (PCR), Sputum Test Code : 8780123 - Laboratory Services | | | |
|---|--|------------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Specimen Type | Container | Volume | Transport |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |

| 3 | Tuberculosis (<i>Mycobacterium tuberculosis</i>) DNA, (PCR), CSF Test Code : 8780125 - Laboratory Services | | | |
|---|---|------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : Tue-Sat 18:00 | | | |
| | Specimen Type | Container | Volume | Transport |
| | CSF | Sterile CSF Tube | 2mL | 2-8°C |

| 4 | Tuberculosis (<i>Mycobacterium tuberculosis</i>) DNA, (PCR), Bronchoalveolar Lavage Test Code : 8780126 - Laboratory Services | | | |
|---|--|---------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Specimen Type | Container | Volume | Transport |
| | Bronchoalveolar Lavage | Biologic Fluid Tube | 5mL | 2-8°C |

| 5 | Tuberculosis (<i>Mycobacterium tuberculosis</i>) DNA, (PCR), Cystic Fluid Test Code : 8780131 - Laboratory Services | | | |
|---|--|---------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Specimen Type | Container | Volume | Transport |
| | Cystic Fluid | Biologic Fluid Tube | 2mL | 2-8°C |

| 6 | Aspergillus DNA (PCR), Sputum Test Code : 8779946 - Laboratory Services | | | |
|---|--|------------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri - Reported : 3 days | | | |
| | Specimen Type | Container | Volume | Transport |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |

| 7 | Aspergillus DNA (PCR), Pus Test Code : 8779950 - Laboratory Services | | | |
|---|---|---------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri - Reported : 3 days | | | |
| | Specimen Type | Container | Volume | Transport |
| | Pus | Biologic Fluid Tube | 2mL | 2-8°C |

| 8 | Candida DNA (PCR) Test Code : 8748101 - Laboratory Services | | | |
|-------------|---|------------------------|--------|----------|
| | Method : MULTIPLEX PCR - Performed : Mon-Fri - Reported : 3 days | | | |
| | Sample Name | Container | Amount | Transfer |
| | Other | Call the Laboratory | | 2-8° C |
| | Aphth Material | Plastic Collection Cup | | 2-8°C |
| | Whole Blood (EDTA) | Lavender (EDTA) | 5mL | 2-8°C |
| | Biopsy Material | Plastic Collection Cup | | 2-8°C |
| Skin Scraps | Plastic Collection Cup | | 2-8°C | |

| 9 | Chlamydia pneumoniae DNA (PCR) Test Code : 8748601 - Laboratory Services | | | |
|------------------------|--|------------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 4 days | | | |
| | Specimen Type | Container | Volume | Transport |
| | Throat Swab | Nucliswab | | 2-8°C |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |
| | Transtracheal Aspirate | Biologic Fluid Tube | 5mL | 2-8°C |
| Bronchoalveolar Lavage | Biologic Fluid Tube | 5mL | 2-8°C | |

| 10 | Legionella pneumophila DNA (PCR) Test Code : 8754201 - Laboratory Services | | | |
|---------------|--|------------------------|--------|-----------|
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Specimen Type | Container | Volume | Transport |
| | Whole Blood (EDTA) (PCR) | Lavender (EDTA) | 2mL | 2-8°C |
| | Bronchoalveolar Lavage | Biologic Fluid Tube | 5mL | 2-8°C |
| | Transtracheal Aspirate | Biologic Fluid Tube | 5mL | 2-8°C |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |
| Tissue Sample | Plastic Collection Cup | | 2-8°C | |

| 11 | Nocardia DNA (PCR) Test Code : 8779914 - Laboratory Services | | | |
|-----|--|------------------------|--------|----------|
| | Method : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Sample Name | Container | Amount | Transfer |
| | Transtracheal Aspirate | Biologic Fluid Tube | 5mL | 2-8°C |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |
| | Pus | Sterile Tube | 2mL | 2-8°C |
| | Tissue Sample | Plastic Collection Cup | | 2-8°C |
| | Bronchoalveolar Lavage | Biologic Fluid Tube | 2mL | 2-8°C |
| CSF | Sterile CSF Tube | 2mL | 2-8°C | |

| 12 | Respiratory Syncytial Virus (RSV) RNA (PCR) Test Code : 8779823 - Laboratory Services | | | |
|------------------------|---|------------------------|--------|-----------|
| | Methodology : REAL TIME PCR - Performed : Mon-Fri 08:00 - Reported : 3 days | | | |
| | Specimen Type | Container | Volume | Transport |
| | Throat Swab | Nucliswab | | 2-8°C |
| | Nasopharyngeal Aspirate | Biologic Fluid Tube | 2mL | 2-8°C |
| | Sputum | Plastic Collection Cup | 2mL | 2-8°C |
| | Nasopharyngeal Swab | Nucliswab | | |
| | Tissue Sample | Plastic Collection Cup | 1gr | 2-8°C |
| | Bronchoalveolar Lavage | Biologic Fluid Tube | 2mL | 2-8°C |
| Transtracheal Aspirate | Biologic Fluid Tube | 2mL | 2-8°C | |

| 13 | Respiratory Panel, PCR Test Code : 8745001 - Laboratory Services | | | |
|------------------------|--|-----------|--------|-----------|
| | Methodology : REAL TIME PCR - Performed : Mon-Fri 08:00 Sat 08:00 - Reported : Tue-Sat 18:00 Mon 18:00 | | | |
| | Specimen Type | Container | Volume | Transport |
| | Throat Swab | Nucliswab | | 2-8°C |
| | Nasopharyngeal Swab | Nucliswab | | 2-8°C |
| Bronchoalveolar Lavage | Biologic Fluid Tube | 5mL | 2-8°C | |

| | | | | |
|-----------|---|------------------------|----------|-----------|
| 14 | Tuberculosis (Mycobacterium tuberculosis) DNA, (PCR), Paraffin Section Test Code : 8780132 - Laboratory Services | | | |
| | Methodology : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week | | | |
| | Specimen Type | Container | Volume | Transport |
| | Paraffin Section | Plastic Collection Cup | 50Mikron | 2-8°C |

| | | | | |
|-----------|---|----------------|--------|----------|
| 15 | Autoimmune Encephalitis Panel, Serum Test Code : 8623530 - Laboratory Services | | | |
| | Method : IFA - Performed : Wed 08:00 - Reported : Fri 18:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | Serum | Separator Tube | 1mL | 2-8°C |

| | | | | |
|-----------|---|------------------|--------|----------|
| 16 | Autoimmune Encephalitis Panel, CSF Test Code : 8623531 - Laboratory Services | | | |
| | Method : IFA - Performed : Wed 08:00 - Reported : Fri 18:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | CSF | Sterile CSF Tube | 1mL | 2-8°C |