

1	Catecholamines (Adr,Norad,Dop.), Urine Random Test Code : 8925694 - Laboratory Services			
	Methodology : LC-MS/MS - Performed : Tue 08:00 - Reported : Wed 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Urine, Random	Plastic Collection Cup	1mL	2-8°C

2	Human Herpes Virus Type 8 (HHV 8) DNA (PCR), Other Material Test Code : 8780194 - Laboratory Services			
	Method : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Other	Call the Laboratory		2-8°C

3	Human Herpes Virus Type 8 (HHV 8) DNA (PCR), Whole Blood Test Code : 8780021 - Laboratory Services			
	Method : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Whole Blood (EDTA) (PCR)	Lavender (EDTA)	2mL	2-8°C

4	Human Herpes Virus Type 8 (HHV 8) DNA (PCR), Serum Test Code : 8780020 - Laboratory Services			
	Method : PCR - Performed : Mon-Fri 08:00 - Reported : 1 week			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	2mL	2-8°C

5	Human Herpes Virus Type 8 (HHV 8) DNA (PCR), Paraffin Section Test Code : 8780019 - Laboratory Services									
	Method : PCR - Performed : Mon 08:00 Fri 08:00 - Reported : 1 week									
	<table border="1"> <thead> <tr> <th>Sample Name</th> <th>Collection Container/Tube</th> <th>Specimen Volume</th> <th>Transfer</th> </tr> </thead> <tbody> <tr> <td>Paraffin Section</td> <td>Plastic Collection Cup</td> <td></td> <td>2-8°C</td> </tr> </tbody> </table>	Sample Name	Collection Container/Tube	Specimen Volume	Transfer	Paraffin Section	Plastic Collection Cup		2-8°C	
Sample Name	Collection Container/Tube	Specimen Volume	Transfer							
Paraffin Section	Plastic Collection Cup		2-8°C							

6	Human Herpes Virus Type 8 (HHV 8) DNA (PCR), Biopsy Material Test Code : 8780018 - Laboratory Services									
	Method : PCR - Performed : Mon 08:00 Fri 08:00 - Reported : 4 days									
	<table border="1"> <thead> <tr> <th>Sample Name</th> <th>Collection Container/Tube</th> <th>Specimen Volume</th> <th>Transfer</th> </tr> </thead> <tbody> <tr> <td>Biopsy Material</td> <td>Plastic Collection Cup</td> <td></td> <td>2-8°C</td> </tr> </tbody> </table>	Sample Name	Collection Container/Tube	Specimen Volume	Transfer	Biopsy Material	Plastic Collection Cup		2-8°C	
Sample Name	Collection Container/Tube	Specimen Volume	Transfer							
Biopsy Material	Plastic Collection Cup		2-8°C							

7	Newborn Metabolic Screening (Comprehensive Panel) Test Code : 8925693 - Laboratory Services									
	Method : - Performed : Mon-Fri 08:00 - Reported : 1 week									
	<table border="1"> <thead> <tr> <th>Sample Name</th> <th>Collection Container/Tube</th> <th>Specimen Volume</th> <th>Transfer</th> </tr> </thead> <tbody> <tr> <td>Capillary Blood</td> <td>Filter Paper</td> <td></td> <td>18-24°C</td> </tr> </tbody> </table>	Sample Name	Collection Container/Tube	Specimen Volume	Transfer	Capillary Blood	Filter Paper		18-24°C	
Sample Name	Collection Container/Tube	Specimen Volume	Transfer							
Capillary Blood	Filter Paper		18-24°C							

Urinary System Stone Disease Panel, Urine (24 h)

Test Code : **8235506** - Laboratory Services

8

Method : **PHOTOMETRY/ENZYMATIC/ISE** - Performed : **Fri 09:00** - Reported : **Tue 18:00**

Sample Name	Collection Container/Tube	Specimen Volume	Transfer
24-hour Urine	Plastic 24-Hour Urine Container	10mL	2-8°C