

Laboratory Test Updates

(effective on December 22, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8616102	C1q Esterase Inhibitor Activity							X			
2	C8666809	Neurofilament Light Chain (NfL)						X				
3	C8406001	Acid Phosphatase			X							
4	C8406601	Acid Phosphatase, Prostatic			X							
5	C8443202	Thyroglobulin (hTg), Washout			X							
6	C8427001	Sex Hormone Binding Globulin (SHBG)			X							
7	C8631604	CA 72-4			X							
8	C8200301	Acetaminophen (Paracetamol)		X	X							
9	C8019601	Salicylate		X	X							
10	C8019801	Theophylline		X	X							
11	C8635301	Interleukin-6 (IL-6)			X							
12	C8266801	Erythropoietin			X							
13	C8631610	Neuron Specific Enolase (NSE), CSF			X							
14	C8631605	Neuron Specific Enolase (NSE)			X							
15	C8217201	Apolipoprotein A1			X							
16	C8217204	Apolipoprotein A2			X							

[illegible]

1	C1q Esterase Inhibitor Activity		
	Test Code : 8616102 - Laboratory Services		
	Methodology : COAGULOMETRY - Performed : Thu 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Plasma (Citrate), Frozen	Blue top (Sodium citrate)	0,5mL
	Transport		<0°C
	Additional Information 		
	The plasma sample should be frozen by portioning it into at least two separate tubes. Unacceptable conditions: Hemolyzed, lipemic, icteric and fibrin.		

2	Neurofilament Light Chain (NfL)		
	Test Code : 8666809 - Laboratory Services		
	Methodology : CLIA - Performed : Wed 12:00 - Reported : 14 days		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1ml
	Transport		2-8°C

3	Acid Phosphatase		
	Test Code : 8406001 - Laboratory Services		
	Methodology : ENZYMATIC - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	0,5mL
	Transport		2-8°C

4	Acid Phosphatase, Prostatic		
	Test Code : 8406601 - Laboratory Services		
	Methodology : ENZYMATIC - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL
	Transport		2-8°C

5	Thyroglobulin (hTg), Washout Test Code : 8443202 - Laboratory Services			
	Methodology : ECLIA - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Biologic Fluid	Biologic Fluid Tube	1mL	2-8°C

6	Sex Hormone Binding Globulin (SHBG) Test Code : 8427001 - Laboratory Services			
	Methodology : ECLIA - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

7	CA 72-4 Test Code : 8631604 - Laboratory Services			
	Methodology : ECLIA - Performed : Mon 08:00 Fri 08:00 - Reported : Tue 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

8	Acetaminophen (Paracetamol) Test Code : 8200301 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Plain Tube	0,3mL	2-8°C

9	Salicylate Test Code : 8019601 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Plain Tube	1mL	2-8°C

10	Theophylline Test Code : 8019801 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Plain Tube	0,5mL	2-8°C

11	Interleukin-6 (IL-6) Test Code : 8635301 - Laboratory Services			
	Methodology : CLIA - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	0,5mL	<0°C

12	Erythropoietin Test Code : 8266801 - Laboratory Services			
	Methodology : CLIA - Performed : Mon 08:00 Fri 08:00 - Reported : Tue 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

13	Neuron Specific Enolase (NSE), CSF Test Code : 8631610 - Laboratory Services			
	Methodology : ECLIA - Performed : Mon 08:00 Fri 08:00 - Reported : Tue 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	CSF	Sterile CSF Tube	1mL	2-8°C

14	Neuron Specific Enolase (NSE) Test Code : 8631605 - Laboratory Services			
	Methodology : ECLIA - Performed : Mon 08:00 Fri 08:00 - Reported : Tue 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	2mL	2-8°C

15	Apolipoprotein A1 Test Code : 8217201 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	0,5mL	2-8°C

16	Apolipoprotein A2 Test Code : 8217204 - Laboratory Services			
	Methodology : NEPHELOMETRY - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

17	Apolipoprotein B Test Code : 8217202 - Laboratory Services			
	Methodology : TURBIDIMETRY - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	0,5mL	2-8°C

18	Fructosamine Test Code : 8298502 - Laboratory Services			
	Methodology : PHOTOMETRY - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

19	Pseudocholinesterase Test Code : 8248001 - Laboratory Services			
	Methodology : ENZYMATIC - Performed : Tue 08:00 Thu 08:00 - Reported : Wed 18:00 Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	0,5mL	2-8°C

