


1

Prolonged Jaundice Panel
 Test Code : **8925704** - Laboratory Services

Method : - Performed : **Mon-Sat 12:00** - Reported : **1 week**

Sample Name	Collection Container/Tube	Specimen Volume	Transfer
Whole Blood (EDTA)	Lavender (EDTA)	2mL	2-8°C
Serum	Separator Tube	2mL	2-8°C
Plasma (NaF)	Gray (Sodium Fluoride)	2mL	2-8°C
Urine, Random	Plastic Collection Cup	1mL	2-8°C
Urine, Sterile	Plastic Collection Cup	2mL	2-8°C
Whole Blood (EDTA)	Lavender (EDTA) (CBC)	2mL	2-8°C

Additional Information 


This panel contains, Whole Blood Count, Thyroid Stimulating Hormone (TSH), Thyroxine, Free (FT4), Urinalysis, Culture, Urine , Bilirubin, Direct, Bilirubin, Total, Glucose-6-Phosphate Dehydrogenase (G6PD), Glucose, Fasting, Pyruvate Kinase. All sample types defined in the test panel are required.

2

Polycystic Ovarian Syndrome (PCOS) Diagnostic Tests
 Test Code : **8818083** - Laboratory Services

Method : - Performed : **Mon 08:00 Thu 08:00** - Reported : **Tue 18:00 Fri 18:00**

Sample Name	Collection Container/Tube	Specimen Volume	Transfer
Serum	Separator Tube	5mL	2-8°C

Additional Information 


Panel Contains; Follicle Stimulating Hormone (FSH), Luteinizing Hormone (LH), 17-Alpha-Hydroxyprogesterone (Direct), Dehydroepiandrosterone Sulphate (DHEA-S), Testosterone, Total, Sex Hormone Binding Globulin (SHBG), Testosterone, Free, Estradiol (E2), Thyroid Stimulating Hormone (TSH), Insulin.

3

Urinary Stone Disease Analysis Panel, Serum
 Test Code : **8235507** - Laboratory Services

Method : **PHOTOMETRY/ISE/CLIA** - Performed : **Mon-Sat 14:00** - Reported : **Mon-Sat 18:00**

Sample Name	Collection Container/Tube	Specimen Volume	Transfer
Serum and Serum, Frozen	Separator Tube	4mL	2-8°C

Additional Information 

Serum and Serum, Frozen
 Test contains; Bicarbonate, Phosphorus, Inorganic, Calcium (Ca), Chloride (Cl), Creatinine, Serum, Potassium (K), Sodium (Na), Urea, Blood Urea Nitrogen (BUN), Uric Acid, Parathyroid Hormone (PTH), Intact. Bicarbonate test; It is recommended to take the blood sample into a gel tube, centrifuge within 2 hours, and transfer it without opening the tube cap after centrifugation.

4

Genital panel (PCR)Test Code : **8749102** - Laboratory ServicesMethod : **MULTIPLEX PCR** - Performed : **Mon 08:00 Wed 08:00** - Reported : **Tue 18:00 Thu 18:00**

Sample Name	Collection Container/Tube	Specimen Volume	Transfer
Urethral Discharge Material	Nucliswab		2-8 °C
Conjunctival Swab	Nucliswab		2-8 °C
Cervical Swab	Nucliswab		2-8 °C
Urine, First Morning Void	Sterile Tube	1,5mL	2-8 °C
Semen	Plastic Collection Cup	1,5mL	2-8 °C
Prostatic fluid	Nucliswab		2-8 °C
Oropharyngeal Swab	Nucliswab		2-8°C
Rectal Swab	Nucliswab		2-8°C

Additional Information

This panel includes; Chlamydia trachomatis, Mycoplasma genitalium, Neisseria gonorrhoeae, Trichomonas vaginalis, Gardnerella vaginalis, Ureaplasma urealyticum, Ureaplasma parvum, Herpes Simplex Virus 1, Herpes Simplex Virus 2, Mycoplasma hominis