

Laboratory Test Updates

(effective on February 09, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8925713	CD4 T-Cell Recent Thymic Emigrants (CD4 RTEs)						X				
2	C8459001	Vitamin A (Retinol)			X							
3	C8444601	Vitamin E (Alpha-Tocopherol)			X							
4	C8029953	PLAC Test (Ip-PLA2)		X			X					
5	C8200902	Acetone/keton, Urine, Qualitative										X
6	C8266801	Erythropoietin			X							

1	<p>CD4 T-Cell Recent Thymic Emigrants (CD4 RTEs)</p> <p>Test Code : 8925713 - Laboratory Services</p>		
	<p>Method : FLOW CYTOMETRY - Performed : Mon-Sat 08:00 - Reported : 3 days</p>		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood	Lavender top (EDTA)	3mL
	<p>Additional Information ▼</p>		
<p>Hemolyzed and lipemic samples are not accepted for analysis.</p>			

2	Vitamin A (Retinol) Test Code : 8459001 - Laboratory Services			
	Method : HPLC - Performed : Wed 08:00 - Reported : Fri 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum, Frozen	Separator Tube	0,5mL	<0°C

3	Vitamin E (Alpha-Tocopherol) Test Code : 8444601 - Laboratory Services			
	Method : HPLC - Performed : Wed 08:00 - Reported : Fri 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum, Frozen	Separator Tube	1mL	<0°C

4	PLAC Test (Ip-PLA2) Test Code : 8029953 - Laboratory Services			
	Methodology : ELISA - Performed : Wed 08:00 - Reported : Thu 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	0,5mL	<0°C

6	Erythropoietin Test Code : 8266801 - Laboratory Services			
	Method : CLIA - Performed : Mon 08:00 Fri 08:00 - Reported : Tue 18:00 Fri 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C