

# Laboratory Test Updates

(effective on March 09, 2026)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8635304	Interleukin-2 (IL-2)					X					
2	C8233001	Calcium (Ca), Ionized				X						
3	C8265704	Arylsulphatase A, Leucocyte				X						
4	C8265703	Beta-Galactosidase, Leucocyte				X						
5	C8308002	Hexosaminidase, Total				X						
6	C8265705	Sphingomyelinase, Leucocyte				X						
7	C8265708	Alpha-Galactosidase, Leucocyte				X						
8	C8296302	Beta-Glucosidase, Leucocyte				X						
9	C8265706	Galactocerebrosidase, Leucocyte				X						
10	C8265711	Alpha-Glucosidase, Leucocyte				X						
11	C8265699	Alpha-Mannosidase, Leukocytes				X						
12	C8265697	Alpha-Iduronidase, Leukocytes				X						
13	C8265698	Alpha-Fucosidase, Leukocytes				X						

<b>1</b>	<b>Interleukin-2 (IL-2)</b>			
	Test Code : 8635304 - Laboratory Services			
	Methodology : <b>FLOW CYTOMETRY</b> - Performed : <b>Tue 08:00</b> - Reported : <b>1 weeks</b>			
	<b>Specimen Type</b>	<b>Collection Container/Tube</b>	<b>Specimen Volume</b>	<b>Transport</b>
	Serum, Frozen	Separator Tube	2mL	<0°C
<b>Reference Ranges</b> <span style="float: right;">▼</span>				
<b>Normal Values</b>				
<b>Gender</b>	<b>Lower Value</b>	<b>Unit</b>		
	<31.25	pg/mL		

<b>2</b>	<b>Calcium (Ca), Ionized</b>			
	Test Code : 8233001 - Laboratory Services			
	Methodology : <b>PHOTOMETRY</b> - Performed : <b>Mon-Sat 12:00</b> - Reported : <b>Mon-Sat 18:00</b>			
	<b>Specimen Type</b>	<b>Collection Container/Tube</b>	<b>Specimen Volume</b>	<b>Transport</b>
	Serum	Separator Tube	1ml	2-8°C
<b>Additional Information</b> <span style="float: right;">▼</span>				
<b>Serum</b>				
The tourniquet should not be applied during blood collection. If necessary, tourniquet application should not exceed 1 minute. To minimise blood contact with air, it is recommended that the sample be collected in a gel tube and transferred after centrifugation without opening the tube cap.				

<b>3</b>	<b>Arylsulphatase A, Leucocyte</b>			
	Test Code : 8265704 - Laboratory Services			
	Methodology : <b>PHOTOMETRY</b> - Performed : <b>Mon 08:00</b> - Reported : <b>Fri 18:00</b>			
	<b>Specimen Type</b>	<b>Collection Container/Tube</b>	<b>Specimen Volume</b>	<b>Transport</b>
	Whole Blood (EDTA)	Lavender top (EDTA)	10mL	2-8°C

4	Beta-Galactosidase, Leucocyte		
	Test Code : 8265703 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
Whole Blood (EDTA)	Lavender top (EDTA)	20mL	Transport
			2-8°C

5	Hexosaminidase, Total		
	Test Code : 8308002 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
Whole Blood (EDTA)	Lavender top (EDTA)	10mL	Transport
			2-8°C

6	Sphingomyelinase, Leucocyte		
	Test Code : 8265705 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
Whole Blood (EDTA)	Lavender top (EDTA)	8mL	Transport
			2-8°C

7	Alpha-Galactosidase, Leucocyte		
	Test Code : 8265708 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
Whole Blood (EDTA)	Lavender top (EDTA)	8mL	Transport
			2-8°C

8	Beta-Glucosidase, Leucocyte		
	Test Code : 8296302 - Laboratory Services		
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00		
	Specimen Type	Collection Container/Tube	Specimen Volume
Whole Blood (EDTA)	Lavender top (EDTA)	8mL	Transport
			2-8°C

9	Galactocerebrosidase, Leucocyte			
	Test Code : 8265706 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (EDTA)	Lavender top (EDTA)	8mL	2-8°C

10	Alpha-Glucosidase, Leucocyte			
	Test Code : 8265711 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (EDTA)	Lavender top (EDTA)	8mL	2-8°C

11	Alpha-Mannosidase, Leukocytes			
	Test Code : 8265699 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (EDTA)	Lavender top (EDTA)	20mL	2-8°C

12	Alpha-Iduronidase, Leukocytes			
	Test Code : 8265697 - Laboratory Services			
	Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (EDTA)	Lavender top (EDTA)	20mL	2-8°C

13

Alpha-Fucosidase, Leukocytes

Test Code : 8265698 - Laboratory Services

Methodology : FLUOROMETRY - Performed : Mon 08:00 - Reported : Fri 18:00

Specimen Type	Collection Container/Tube	Specimen Volume	Transport
Whole Blood (EDTA)	Lavender top (EDTA)	20mL	2-8°C