

# Laboratory Test Updates

(effective on June 02, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8215008	Dickkopf 3 (DKK3)						X				
2	C8925661	Acanthamoeba DNA (PCR), Biopsy Material										X
3	C8925660	Acanthamoeba DNA (PCR), CSF										X
4	C8818077	Daratumumab Interference Test (MALDI-TOF)									X	
5	C8623526	CIDP / Nodo Paranodopathy						X				

1	<b>Dickkopf 3 (DKK3)</b> Test Code : 8215008 - Laboratory Services			
	Method : ELISA - Performed : Mon 08:00 Wed 08:00 - Reported : 10 days			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Urine, Frozen	Plastic Collection Cup	2mL	<0°C
	Additional Information <span style="float: right;">▼</span> It is suitable for analysis in spot urine samples. The specimen should be frozen and transported in a frozen state. If possible, the first morning urine is recommended, as it is more concentrated and thus increases analytical sensitivity.			

4	<b>Daratumumab Interference Test (MALDI-TOF)</b> Test Code : 8818077 - Laboratory Services			
	Method : NMR - Performed : Wed 08:00 - Reported : 6 days			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	5mL	2-8°C

5	<b>CIDP / Nodo Paranodopathy</b> Test Code : 8623526 - Laboratory Services			
	Methodology : IFA - Performed : Wed 08:00 - Reported : 2 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	0,5mL	2-8°C
	Plasma (EDTA)	Lavender (EDTA)	0,5mL	2-8°C
	CSF	Sterile CSF Tube	0,5mL	2-8°C
Clinical Information <span style="float: right;">▼</span>				
Used in the diagnosis of polyneuropathy.				