

# Laboratory Test Updates

(effective on November 27, 2024)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8029902	Methotrexate							X			
2	C8636001	CD4/CD8 Ratio							X			

<b>1</b>	<b>Methotrexate</b> Test Code : <b>8029902</b> - Laboratory Services			
	Methodology : <b>EMIT</b> - Performed : <b>Mon-Sat</b> - Reported : <b>4 hours</b>			
	Specimen Type	Container	Volume	Transport
	Serum	Plain Tube	1mL	2-8°C
	Additional Information <span style="float: right;">▼</span>			
Information about the blood collecting time after the administration of medication should be added. (etc.24h 48 h 96 h)				
<b>Serum</b> Unacceptable Conditions: Separator tubes. Protect from light.				

2

### CD4/CD8 Ratio

Test Code : **8636001** - Laboratory Services

Methodology : **FLOW CYTOMETRY** - Performed : **Mon 08:00 Tue 08:00 Wed 08:00 Thu 08:00 Fri 08:00 Sat 08:00** -  
Reported : **Wed 18:00 Thu 18:00 Fri 18:00 Mon 18:00 Tue 18:00 Tue 18:00**

Specimen Type	Container	Volume	Transport
Whole Blood	Lavender (EDTA)	5mL	18-24°C
Biologic Fluid	Lavender (EDTA)	6mL	18-24°C
CSF	Lavender (EDTA)	6mL	18-24°C

#### Additional Information

Please contact with the laboratory before collecting the sample.Fresh sample is required.It is recommended that specimens arrive within 24 hours of draw.Test contains CD45, CD3, CD4, CD8.

