

Laboratory Test Updates

(effective on December 29, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8017001	Gentamycin				X						
2	C8305001	Methemoglobin							X			

1	Gentamycin		
	Test Code : 8017001 - Laboratory Services		
	Methodology : IMMUNOASSAY - Performed : Mon 12:00 Wed 12:00 - Reported : 8 days		
	Specimen Type	Collection Container/Tube	Specimen Volume
	Serum	Separator Tube	1mL
	Serum	Plain Tube	1mL
			Transport
			2-8°C
			2-8°C

2	Methemoglobin		
	Test Code : 8305001 - Laboratory Services		
	Method : SPECTROPHOTOMETRY - Performed : Mon-Fri 08:00 - Reported : Tue-Sat 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (EDTA)	Lavender top (EDTA)	4mL
			Transfer
			2-8°C
	Additional Information		▼
	<p>The sample taken for the methemoglobin test must be delivered to the Central Laboratory within 72 hours at 2-8°C. Samples that have been stored at room temperature will be rejected.</p>		

