

| | | | | |
|---|--|------------------|--------|----------|
| 1 | Culture, CSF Test Code : 8707507 - Laboratory Services | | | |
| | Method : C/I/S - Performed : Mon-Sat 12:00 - Reported : 6 days 17:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | CSF | Sterile CSF Tube | 0,3mL | 18-24°C |

| | | | | |
|-----------------|---|------------------------|--------|----------|
| 2 | Chlamydia trachomatis Antigen Test Code : 8727001 - Laboratory Services | | | |
| | Method : IMMUNOASSAY - Performed : Mon-Fri 08:00 Sat 08:00 - Reported : Tue-Sat 18:00 Mon 18:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | Conjunctival Swab | Dacron/ Rayon Swab | | 2-8°C |
| | Urine | Plastic Collection Cup | 20mL | 2-8°C |
| | Cervical Swab | Dacron/ Rayon Swab | | 2-8°C |
| | Urethral Discharge Material | Dacron/ Rayon Swab | | 2-8°C |
| Prostatic fluid | Dacron/ Rayon Swab | | 2-8°C | |

| | | | | |
|--------|---|--------------------|--------|----------|
| 3 | Influenza Antigen (A and B) Test Code : 8740001 - Laboratory Services | | | |
| | Method : IMMUNOASSAY - Performed : Mon-Sat - Reported : 2 hours | | | |
| | Sample Name | Container | Amount | Transfer |
| | Throat Swab | Dacron/ Rayon Swab | | 2-8°C |
| | Nasal Discharge | Dacron/ Rayon Swab | | 2-8°C |
| | Nasopharyngeal Aspirate | Dacron/ Rayon Swab | | 2-8°C |
| | Nasopharyngeal Swab | Dacron/ Rayon Swab | | 2-8°C |
| Sputum | Plastic Collection Cup | | 2-8°C | |

| | | | | |
|---|---|------------------------|--------|----------|
| 4 | Culture, Throat Test Code : 8706001 - Laboratory Services | | | |
| | Method : C/I - Performed : Mon-Sat 12:00 - Reported : 2 days 17:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | Throat Swab | Culture Transport Swab | | 18-24°C |

| 5 | Gonococcus Culture Test Code : 8708101 - Laboratory Services | | | |
|-----------------------------|--|------------------------|---------|----------|
| | Method : <i>C/VS</i> - Performed : Mon-Sat 12:00 - Reported : 4 days 17:00 | | | |
| | Sample Name | Container | Amount | Transfer |
| | Vaginal Discharge | Culture Transport Swab | | 18-24°C |
| | Cervical Swab | Culture Transport Swab | | 18-24°C |
| Urethral Discharge Material | Culture Transport Swab | | 18-24°C | |

| 6 | ESBL Positive/ Carbapenem Resistant Organism Screening Test Code : 8708147 - Laboratory Services | | | |
|--------------------|--|------------------------|---------|-----------|
| | Methodology : <i>C/VS</i> - Performed : Mon-Fri 12:00 - Reported : 3 days 17:00 | | | |
| | Specimen Type | Container | Volume | Transport |
| | Rectal Swab | Culture Transport Swab | | 18-24°C |
| | Other | Call the Laboratory | | 18-24°C |
| Oropharyngeal Swab | Culture Transport Swab | | 18-24°C | |

| 7 | Culture, Legionella Test Code : 8708102 - Laboratory Services | | | |
|--------|--|------------------------|--------|-----------|
| | Methodology : CULTURE - Performed : Mon-Sat 12:00 - Reported : 7 days 17:00 | | | |
| | Specimen Type | Container | Volume | Transport |
| | Urine, Sterile | Plastic Collection Cup | 5mL | 2-8°C |
| | Nasopharyngeal Swab | Culture Transport Swab | | 2-8°C |
| | Bronchoalveolar Lavage | Sterile Tube | 2mL | 2-8°C |
| | Transtracheal Aspirate | Sterile Tube | 2mL | 2-8°C |
| | Tissue Sample | Sterile Tube | | 2-8°C |
| | Biologic Fluid | Sterile Tube | 2mL | 2-8°C |
| Sputum | Sterile Tube | | 2-8°C | |

Culture, FungiTest Code : **8710201** - Laboratory Services**8**Methodology : **C/I/S** - Performed : **Mon-Sat 12:00** - Reported : **21 days 17:00**

| Specimen Type | Container | Volume | Transport |
|------------------------|------------------------|--------|-----------|
| Nail Scraps | Plastic Collection Cup | | 18-24°C |
| Urine, Sterile | Plastic Collection Cup | | 2-8°C |
| Throat Swab | Culture Transport Swab | | 2-8°C |
| Skin Scraps | Plastic Collection Cup | | 18-24°C |
| Vaginal Discharge | Culture Transport Swab | | 2-8°C |
| Hair | Plastic Collection Cup | | 18-24°C |
| Other | Call the Laboratory | | 2-8°C |
| Whole Blood | Blood Culture Bottle | 10mL | 18-24°C |
| Sputum | Plastic Collection Cup | | 2-8°C |
| Peripheral Stem Cell | Sterile Glass Tube | 1mL | 18-24°C |
| Bronchoalveolar Lavage | Biologic Fluid Tube | 1mL | 2-8°C |

9

Culture, Aerobic (All Sample Types)

Test Code : 8707501 - Laboratory Services

Methodology : CA/S - Performed : Mon-Sat 12:00 - Reported : 6 days 17:00

| Specimen Type | Container | Volume | Transport |
|-----------------------------|------------------------|--------|-----------|
| Skin Material | Call the Laboratory | | 2-8°C |
| Surgical Material | Call the Laboratory | | 2-8°C |
| Pleural Fluid | Biologic Fluid Tube | 1mL | 2-8°C |
| Biologic Fluid | Biologic Fluid Tube | 1mL | 2-8°C |
| Nail Scraps | Plastic Collection Cup | | 2-8°C |
| Cervical Swab | Biologic Fluid Tube | | 2-8°C |
| Cord Blood | Plain Tube | | 2-8°C |
| Pericardial Fluid | Biologic Fluid Tube | 1mL | 2-8°C |
| Ear Discharge | Biologic Fluid Tube | | 2-8°C |
| Tissue Sample | Plastic Collection Cup | | 2-8°C |
| Rectal Swab | Culture Transport Swab | | 2-8°C |
| Nasopharyngeal Swab | Culture Transport Swab | | 2-8°C |
| Nasal Discharge | Culture Transport Swab | | 2-8°C |
| Urethral Discharge Material | Culture Transport Swab | | 2-8°C |
| Peritoneal Fluid | Biologic Fluid Tube | 1mL | 2-8°C |
| Cord | Plastic Collection Cup | | 2-8°C |
| Gastric Juice | Biologic Fluid Tube | 1mL | 2-8°C |
| Ophthalmic Discharge | Biologic Fluid Tube | 0,5mL | 2-8°C |
| Biopsy Material | Plastic Collection Cup | | 2-8°C |
| Joint Fluid | Biologic Fluid Tube | 1mL | 2-8°C |
| Other | Call the Laboratory | | 2-8°C |
| Cystic Fluid | Sterile Tube | | 2-8°C |
| Bone Marrow | Biologic Fluid Tube | | 2-8°C |
| Vulva | Culture Transport Swab | | 2-8°C |
| Prostatic fluid | Biologic Fluid Tube | | 2-8°C |

| | | | | |
|---|--|------------------------|---------------|------------------|
| 10 | Occult Blood, Stool Test Code : 8227001 - Laboratory Services | | | |
| | Methodology : IMMUNOCHROMATOGRAPHY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00 | | | |
| | Specimen Type | Container | Volume | Transport |
| | Stool | Plastic Collection Cup | 1gr | 2-8°C |
| | Additional Information ▼ | | | |
| <p>Stool</p> <p>Except for compulsory cases, aspirin or other anti-rheumatic drugs, painkillers, antipyretics should not be taken 1 week before the test, or the physician should be informed of this situation. Suppository drugs should not be used 2 days before the test. If the person who will have the test is in her menstrual period, there will be blood in her urine. If the patient has a bleeding hemorrhoid or anal fissure and has difficulty defecating, the patient should inform the laboratory and, if possible, have the test done after these conditions have passed. Samples can be stored at room temperature for 2 hours or at 2-8°C for 3 days.</p> | | | | |

| | | | | |
|-----------|--|------------------|---------------|------------------|
| 11 | Galactose, Free Test Code : 8276001 - Laboratory Services | | | |
| | Methodology : PHOTOMETRY - Performed : Tue 08:00 - Reported : 2 weeks | | | |
| | Specimen Type | Container | Volume | Transport |
| | Capillary Blood | Filter Paper | | 18-24°C |