

(effective on November 3, 2025)

[illegible]

1	Copper [Mass/time] in 24 hour Urine			
	Test Code : 8252506 - Laboratory Services			
	Method : AAS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	24-hour Urine, HCl treated	24-Hour Urine Container	10mL	2-8°C

2	Copper/Creatinine [Mass Ratio] in Urine			
	Test Code : 8252505 - Laboratory Services			
	Method : AAS - Performed : Tue 08:00 Wed 08:00 Fri 08:00 - Reported : Thu 18:00 Fri 18:00 Tue 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Urine, Random	Plastic Collection Cup	10mL	2-8°C

3	Bleeding Time (PFA)			
	Test Code : 8500202 - Laboratory Services			
	Methodology : AGREGOMETRY - Performed : Mon-Sat 12:00 - Reported : Mon-Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Whole Blood (Citrate)	Blue top (Sodium citrate)	2mL	18-24°C
	Additional Information 			
	This test is available at Acibadem Altunizade, Maslak, and Atakent Hospitals. A fresh specimen is required for analysis. Due to specimen stability, the sample should be collected close to the courier dispatch time and must reach the laboratory within 4 hours of collection.			