

Laboratory Test Updates

(effective on November 17, 2025)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8666803	Francisella tularensis (Tularemia) Antibody, IgM						X				
2	C8684906	Schistosoma Antibody		X		X					X	
3	C8674701	Parvovirus B19 Antibody, IgG		X								
4	C8674702	Parvovirus B19 Antibody, IgM		X								
5	C8660906	Borrelia burgdorferi Antibody, IgG (Lyme Disease)		X								
6	C8661802	Borrelia burgdorferi Antibody, IgM (Lyme Disease)		X								
7	C8662205	Brucella Antibody, IgG		X								
8	C8662204	Brucella Antibody, IgM		X								
9	C8019703	ALEX3, Multiplex Comprehensive Allergy Panel (300 Parameters)	X									
10	C8749108	Enteric Panel			X							
11	C8352022	TNF-Alpha		X			X					
12	C8749157	Phenome Microbiota (Metagenom)	X									
13	C8208804	Aldosterone/ Plasma Renin Activity (PRA) Ratio				X			X			
14	C8233001	Calcium (Ca), Ionized										X
15	C8424401	Renin, Direct (Active Renin)					X					
16	C8557602	Platelet Function Assay Test, LTA									X	
17	C8684910	Interferon gamma		X		X						

1	Francisella tularensis (Tularemia) Antibody, IgM Test Code : 8666803 - Laboratory Services			
	Method : ELISA - Performed : Tue 08:00 Thu 08:00 - Reported : 7 days			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

2	Schistosoma Antibody Test Code : 8684906 - Laboratory Services			
	Methodology : IHA/ IFT - Performed : Tue 08:00 - Reported : 2 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	1mL	<0°C

3	Parvovirus B19 Antibody, IgG Test Code : 8674701 - Laboratory Services			
	Method : EIA, CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

4	Parvovirus B19 Antibody, IgM Test Code : 8674702 - Laboratory Services			
	Methodology : EIA, CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

5	Borrelia burgdorferi Antibody, IgG (Lyme Disease) Test Code : 8660906 - Laboratory Services			
	Methodology : EIA, CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

6	Borrelia burgdorferi Antibody, IgM (Lyme Disease) Test Code : 8661802 - Laboratory Services			
	Method : EIA, CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

7	Brucella Antibody, IgG Test Code : 8662205 - Laboratory Services			
	Method : EIA, CLIA - Performed : Tue 08:00 Fri 08:00 - Reported : Wed 18:00 Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	1mL	2-8°C

8	Brucella Antibody, IgM Test Code : 8662204 - Laboratory Services			
	Methodology : EIA, CLIA - Performed : Tue 08:00 Fri 08:00 - Reported : Wed 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum	Separator Tube	1mL	2-8°C

9	ALEX3, Multiplex Comprehensive Allergy Panel (300 Parameters) Test Code : 8019703 - Laboratory Services			
	Method : EIA - Performed : Fri 08:00 - Reported : Mon 18:00			
	Sample Name	Collection Container/Tube	Specimen Volume	Transfer
	Serum	Separator Tube	0,5mL	2-8°C

10	Enteric Panel Test Code : 8749108 - Laboratory Services			
	Methodology : REAL TIME PCR - Performed : Mon-Fri 08:00 Sat 08:00 - Reported : Tue-Sat 18:00 Sat 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Stool	Plastic Collection Cup	1gr	2-8°C

11	TNF-Alpha Test Code : 8352022 - Laboratory Services			
	Methodology : CLIA - Performed : Tue 08:00 - Reported : 1 week			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Serum, Frozen	Separator Tube	1mL	<0°C
	Reference Ranges 			
	Normal Values			
	Gender	Lower Value	Unit	
		<13.7	pg/mL	

12	Phenome Microbiota (Metagenom) Test Code : 8749157 - Laboratory Services			
	Methodology : NEXT GENERATION SEQUENCING (NGS) - Performed : Mon-Fri - Reported : 5-6 weeks			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Apply to Laboratory	Special Sample Collection Kit		18-24°C

13	Aldosterone/ Plasma Renin Activity (PRA) Ratio			
	Test Code : 8208804 - Laboratory Services			
	Methodology : CLIA - Performed : Tue 12:00 - Reported : Fri 12:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (EDTA), Frozen	Lavender top (EDTA)	2mL	<0°C

Additional Information 

Plasma (EDTA), Frozen
The plasma sample should be frozen by portioning it into at least two separate tubes.

15	Renin, Direct (Active Renin)			
	Test Code : 8424401 - Laboratory Services			
	Methodology : CLIA - Performed : Mon 08:00 Wed 08:00 Fri 08:00 - Reported : Tue 18:00 Thu 18:00 Mon 18:00			
	Specimen Type	Collection Container/Tube	Specimen Volume	Transport
	Plasma (EDTA), Frozen	Lavender top (EDTA)	1mL	<0°C

Reference Ranges 

Special Values				
Gender	Description	Lower Value	Upper Value	Unit
	Upright	5.3	99.1	uIU/mL
	Supine	4.2	59.7	uIU/mL

16	Platelet Function Assay Test, LTA		
	Test Code : 8557602 - Laboratory Services		
	Method : LTA - Performed : Mon-Fri 16:00 Sat 12:00 - Reported : Tue-Sat 18:00 Mon 18:00		
	Sample Name	Collection Container/Tube	Specimen Volume
	Whole Blood (Citrate)	Blue top (Sodium citrate)	6mL
			Transfer
			18-24°C

17

Interferon gamma			
Test Code : 8684910 - Laboratory Services			
Method : ELISA - Performed : Tue 08:00 - Reported : 1 week			
Sample Name	Collection Container/Tube	Specimen Volume	Transfer
Serum, Frozen	Separator Tube	2mL	<0°C
Reference Ranges			
Normal Values			
Gender	Lower Value	Unit	
	<0.9	IU/mL	