

Laboratory Test Updates

(effective on August 05, 2024)

	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8818082	Mucopolysaccharides, TLC, Urine (random)						X				
2	C8779905	HPV DNA Typing (PCR)										X
3	C8779906	HPV DNA, Screening Test (PCR)										X
4	C8780008	HPV DNA Typing (PCR), Biopsy Material									X	
5	C8780009	HPV DNA Typing (PCR), Urine									X	
6	C8780010	HPV DNA Typing (PCR), Liqui-Prep									X	
7	C8780011	HPV DNA Typing (PCR), Paraffin Section									X	
8	C8780012	HPV DNA Typing (PCR), Swab									X	
9	C8780013	HPV DNA Typing (PCR), Thin Prep									X	
10	C8780004	HPV DNA, Screening Test (PCR), Biopsy Material									X	
11	C8780005	HPV DNA, Screening Test (PCR), Urine									X	
12	C8780006	HPV DNA, Screening Test (PCR), Swab									X	
13	C8780007	HPV DNA, Screening Test (PCR), Thin Prep									X	
14	C8707502	Culture, Pus			X							
15	C8707501	Culture, Aerobic (All Sample Types)			X							
16	C8707511	Culture, Anaerobic			X							

1	Mucopolysaccharides, TLC, Urine (random) Test Code : 8818082 - Laboratory Services			
	Method : CHROMATOGRAPHY - Performed : Thu 08:00 - Reported : 2 weeks			
	Sample Name	Container	Amount	Transfer
	Urine, Random, Frozen	Plastic Collection Cup	5mL	<0°C

4	HPV DNA Typing (PCR), Biopsy Material Test Code : 8780008 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 12:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Biopsy Material	Plastic Collection Cup		2-8°C
Additional Information ▼				
This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.				


5	HPV DNA Typing (PCR), Urine Test Code : 8780009 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 12:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Urine	Sterile Tube	5mL	2-8°C
Additional Information ▼				
This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.				

6

HPV DNA Typing (PCR), Liqui-Prep
 Test Code : **8780010** - Laboratory Services

Method : **PCR, DNA SEQUENCING** - Performed : **Mon-Fri 12:00** - Reported : **10 days**

Sample Name	Container	Amount	Transfer
Liqui-Prep	Call the Laboratory		2-8°C

Additional Information 

This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.

7

HPV DNA Typing (PCR), Paraffin Section
 Test Code : **8780011** - Laboratory Services

Method : **PCR, DNA SEQUENCING** - Performed : **Mon-Fri 12:00** - Reported : **10 days**

Sample Name	Container	Amount	Transfer
Paraffin Section	Plastic Collection Cup		2-8°C

Additional Information 


This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.

8

HPV DNA Typing (PCR), Swab
 Test Code : **8780012** - Laboratory Services

Method : **PCR, DNA SEQUENCING** - Performed : **Mon-Fri 12:00** - Reported : **10 days**

Sample Name	Container	Amount	Transfer
Swab	Nucliswab		2-8°C

Additional Information 

This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.

Swab
 HPV samples should be taken as skin swab/ scrapings and transferred in nucliswab media.

9	HPV DNA Typing (PCR), Thin Prep Test Code : 8780013 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 12:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Thin Prep	ThinPrep		2-8°C
	Additional Information ▼			
This test includes; HPV 16, HPV 18, HPV 31, HPV 33, HPV 35, HPV 39, HPV 45, HPV 51, HPV 52, HPV 56, HPV 58, HPV 59, HPV 66 ve HPV 68 types.				

10	HPV DNA, Screening Test (PCR), Biopsy Material Test Code : 8780004 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 08:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Biopsy Material	Plastic Collection Cup		2-8°C

11	HPV DNA, Screening Test (PCR), Urine Test Code : 8780005 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 08:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Urine	Sterile Tube	2mL	2-8°C

12	HPV DNA, Screening Test (PCR), Swab Test Code : 8780006 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 08:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Swab	Swap		2-8°C

13	HPV DNA, Screening Test (PCR), Thin Prep Test Code : 8780007 - Laboratory Services			
	Method : PCR, DNA SEQUENCING - Performed : Mon-Fri 08:00 - Reported : 10 days			
	Sample Name	Container	Amount	Transfer
	Thin Prep	ThinPrep		2-8°C

14	Culture, Pus Test Code : 8707502 - Laboratory Services			
	Method : C/I/S - Performed : Mon-Sat 12:00 - Reported : 6 days 17:00			
	Sample Name	Container	Amount	Transfer
	Pus	Biologic Fluid Tube	2gr	2-8°C

15

Culture, Aerobic (All Sample Types)

Test Code : 8707501 - Laboratory Services

Methodology : C/I/S - Performed : Mon-Sat 12:00 - Reported : 6 days 17:00

Specimen Type	Container	Volume	Transport
Skin Material	Call the Laboratory		2-8°C
Surgical Material	Call the Laboratory		2-8°C
Pleural Fluid	Biologic Fluid Tube	1mL	2-8°C
Biologic Fluid	Biologic Fluid Tube	1mL	2-8°C
Nail Scraps	Plastic Collection Cup		2-8°C
Cervical Swab	Biologic Fluid Tube		2-8°C
Cord Blood	Plain Tube		2-8°C
Pericardial Fluid	Biologic Fluid Tube	1mL	2-8°C
Ear Discharge	Biologic Fluid Tube		2-8°C
Tissue Sample	Plastic Collection Cup		2-8°C
Rectal Swab	Swap		2-8°C
Nasopharyngeal Swab	Swap		2-8°C
Nasal Discharge	Swap		2-8°C
Urethral Discharge Material	Swap		2-8°C
Peritoneal Fluid	Biologic Fluid Tube	1mL	2-8°C
Cord	Plastic Collection Cup		2-8°C
Gastric Juice	Biologic Fluid Tube	1mL	2-8°C
Ophthalmic Discharge	Biologic Fluid Tube	0,5mL	2-8°C
Biopsy Material	Plastic Collection Cup		2-8°C
Joint Fluid	Biologic Fluid Tube	1mL	2-8°C
Other	Call the Laboratory		2-8°C
Cystic Fluid	Sterile Tube		2-8°C
Bone Marrow	Biologic Fluid Tube		2-8°C
Vulva	Swap		2-8°C
Prostatic fluid	Biologic Fluid Tube		2-8°C

16

Culture, Anaerobic

Test Code : 8707511 - Laboratory Services

Methodology : C/I/S - Performed : Mon-Sat 12:00 - Reported : 6 days 17:00

Specimen Type	Container	Volume	Transport
Bronchoalveolar Lavage	Biologic Fluid Tube	5mL	18-24°C
Other	Call the Laboratory	5mL	18-24°C
Transtracheal Aspirate	Biologic Fluid Tube	5mL	18-24°C
Biologic Fluid	Biologic Fluid Tube	5mL	18-24°C
Pus	Biologic Fluid Tube	5mL	18-24°C