

Laboratory Test Updates

(effective on May 7, 2024)


	Test Code	Test Name	Name Change	Methodology	Performed/Reported Schedule	Specimen Requirements	Reference Interval	New Tests	Additional Information	Internalized Tests	Reactivated Tests	Inactive
1	C8029953	PLAC test (Ip-PLA2)			X						X	
2	C8661112	Pediatri Food Mix Allergy Panel						X				

1

PLAC test (Ip-PLA2)
 Test Code : **8029953** - Laboratory Services

Method : **ENZYMATIC** - Performed : **Tue 08:00** - Reported : **Wed 18:00**

Sample Name	Container	Amount	Transfer
Serum, Frozen	Separator Tube	0,5mL	<0°C

Additional Information 

Unacceptable Conditions: Grossly hemolyzed, icteric and lipemic specimens.

Serum, Frozen
 Hemolysed sample is not acceptable.

2

Pediatri Food Mix Allergy Panel
 Test Code : **8661112** - Laboratory Services

Methodology : - Performed : **Mon-Thu 08:00 Fri 08:00** - Reported : **Tue-Fri 18:00 Mon 18:00**

Specimen Type	Container	Volume	Transport
Serum	Separator Tube	3mL	2-8°C

Additional Information 

Test contains; Specific IgE- Childrens Food Panel (fx5), Specific IgE- Food Mix 1 (fx1), Specific IgE- Food mix (fx13), Specific IgE- Food Mix 2 (fx2), Specific IgE- Food Mix 3 (fx3), Specific IgE- Food mix (fx50), Specific IgE- Food Mix 6 (fx6), Specific IgE- Food mix (fx73), Spesific IgE-Food Mix (fx15).